

PACKING LIST

ForensiKit: Printing with Magnetic Powder

- One magnetic applicator wand**
This tool is used to hold the magnetic powder as powder is applied to the surface
- One 1-ounce jar magnetic powder**
Black magnetic powder reveals latent prints on light surfaces
- One pack of 2 x 6 inch lift tape**
Use lifting tape to transfer the latent print from the surface to the latent print card
- 25 latent fingerprint cards**
Place your lifted print on the glossy card back.
- One dust mask**
Protect yourself from particulates
- One pair of Nitrile gloves**
Protect against cross-contamination and stains
- Four sample practice surfaces**
 - 1 glass slide
 - 1 sample of unpainted wood
 - 1 sample of primed wood
 - 1 sample of paper

Not included

- Plastic bag or tarp to protect the work area

ForensiKit

Lorem ipsum

Tips & Safety

Printing with Magnetic Powder

- Stop dusting as soon as you see a latent print to avoid removing too much of the magnetic powder.
 - If a very greasy print shows up as a solid blob, additional sweeps with the wand may help.
 - For best results, always use the magnetic applicator wand with the magnet end down.
 - To create fingerprints, rub the side of your nose, forehead or hair, and then touch the surface where you want to leave the print.
-

Work Station Best Practices

- Use on a surface that can be cleaned with water or discarded.
 - Use outdoors if possible.
 - Cover your work area with an empty plastic trash bag, and discard everything after use.
 - Do not use near food or food preparation areas.
 - When finished, clean with copious water.
 - Use protective gear: gloves, breathing mask.
 - Store your forensic supplies in a cool, dry place.
-