

Printing with Fluorescent Magnetic Powder

PACKING LIST

- One magnetic applicator wand
- One 1-oz vial of fluorescent magnetic powder
- One 400 nm UV light
- One pair of orange UV safety glasses
- One dust mask
- One pair of Nitrile gloves

ForensiKit

Tips & Safety

- Choose a work surface that can be cleaned with water or discarded.
 - » Work outdoors if possible.
- Protect your work area with an empty plastic trash bag, and discard everything after use.
- Do not use fingerprint powder near food or food preparation areas.
- Always wear a dust mask and gloves when working with fingerprint powder.
- Store your forensic supplies in a cool, dry place.

Working with magnetic powder

- Stop dusting as soon as you see a latent print to avoid removing too much of the magnetic powder.
- If a very greasy print shows up as a solid blob, additional sweeps with the wand may help.
- For best results, always use the magnetic applicators with the magnet end down.
- If you want to create a print and your fingers are dry, rub the side of your nose, your forehead or your hair to pick up additional natural oils.
- When creating a fingerprint, don't press down. Simply touch the surface where you want to leave the print