

Gunshot Evidence

PACKING LIST

- Four Forenstix-Lead tests
- Three spent cartridges
- One 5X magnifier linen tester
- One 6-inch wood ruler
- One caliper
- One pair of nitrile gloves
- One pair of cotton gloves

ForensiKit

Tips & Safety

The shell casings and Forenstix test paper contain traces of lead and other chemicals. This warrants a handling warning as follows.

- » **WARNING:** This product contains a chemical known to the State of California to cause cancer.
- » **WARNING:** This product can expose you to chemicals including arsenic, which is known to the State of California to cause cancer. For more information, go to www.P65Warnings.ca.gov.

When you perform the Forenstix test for lead be sure to hold the cap closed while squeezing the vial to prevent the materials from escaping. Operate away from your face.

You can test other surfaces with your Forenstix lead tests. Moisten the swab and rub on the suspect area.

IMPORTANT

There are trace amounts of lead in the shell casings. Handle shell casings with gloves while wearing a mask. Avoid contact.

For more information, activities and video demos, visit <https://forensikit.com/gunshot-evidence/>